


# ACTION<sup>®</sup>

*Making a Material Difference<sup>®</sup>*




**Wheelchair Cushions  
Mattress Overlays • Accessories**

# Making a *Material* Difference®

Leading the industry in advanced polymer technology for over 35 years, Action® products offer simple, effective solutions for your skin protection challenges. Our proprietary AKTON® polymer relieves pressure, reduces shear forces, manages heat & humidity build-up and is 100% hassle-free. Its unique chemical structure means you don't have to worry about leaks, cuts, punctures, oozing, deflating or bottoming out. Enjoy peace of mind knowing this ultrasoft polymer provides comfortable, full-memory support that's proven clinically effective and long-lasting durability that is simple and affordable.


## Basic Comfort and Skin Protection

### Action® Pilot™




- 1" Akton® polymer
- moderate risk user
- low profile design
- sizes 16" x 16" to 20" x 20"
- custom sizes available
- 5 year warranty

### Action® Centurian™




- 5/8" foam base w. 5/8" Akton® polymer topper
- low to moderate risk user
- low profile design
- sizes 16" x 16" to 20" x 20"
- custom sizes available to 24" max
- 5 year warranty

### Action® Professional™



- 1 5/8" Akton® polymer
- moderate to high risk user
- low profile design
- sizes 16" x 16" to 20" x 20"
- custom sizes available
- 5 year warranty

### Action® Commuter™



- 1 3/4" foam base w. 1/2" Akton® polymer topper
- low to moderate risk user
- lightweight, portable design
- sizes 16" x 16" to 20" x 20"
- custom sizes available to 24" max
- 5 year warranty

**Cover Options:** All above products sold standard with Basic cover. Incontinent and replacement covers also available as defined in price sheet.

## Action® Mattress Overlays

### HCPCS E0185, Gel-Like Mattress



- 100% Akton® polymer
- moderate to high risk user
- hassle free care & use
- bed hitch option available
- multiple sizes available in 5/8" or 7/8" thickness
- 5 year warranty



800.228.7763 • www.actionproducts.com

Action Products, Inc. • 954 Sweeney Drive •



CE

MADE IN USA





Hassle-Free • Clinically Proven • Repairable • Latex Free


## Advanced Seating & Positioning

### Xact® Classic




- contoured foam base w. visco foam ischial insert
- 9/16" cube Akton® polymer topper
- moderate to high risk user
- sizes 10" x 10" to 20" x 20"
- custom sizes available to 24" max
- 2 year warranty
- HCPCS E2607, Skin Protection & Positioning

### Xact® Lite




- 2 3/4" foam base w. 3/16" pebbled Akton® polymer topper
- active, mobile user
- extreme lightweight, low profile design
- sizes 16" x 16" to 20" x 20"
- 2 year warranty
- HCPCS E2603, Skin Protection

### Xact® Individual



- fully customized foam/ Akton® polymer cushion
- thousands of configurations possible
- therapist input recommended for proper fit to user
- 1/2", 1", 1 1/2" or adjustable angle drop seat option available
- 2 year warranty
- HCPCS E2609, Custom Fabricated

### Xact® Soft



- 3 1/2" contoured foam base
- 9/16" cube Akton® polymer ischial insert
- moderate to high risk user
- soft, maximum immersion design
- sizes 16" x 16" to 20" x 20"
- 2 year warranty
- HCPCS E2607, Skin Protection & Positioning

## Adaptive Use

**Product Features:** • 100% Akton® polymer designs • easily cut & resealed for spot protection • provides comfort and can enhance other cushions  
 • perfect for custom configurations • not intended for use as cushion for at-risk user • repair/reseal kits available



### Cube Surface

- 9/16" cubed Akton® polymer construct
- reduces shear
- contours easily
- sizes 16" x 16" to 20" x 20"
- custom sizes available to 24"


### Flat Surface

- 1/4", 1/2" or 5/8" flat Akton® polymer
- easily cut for custom/ spot protection
- size range 16" x 16" to 20" x 20" \*
- custom sizes available to 24"

### Pebbled Surface

- 1/4", 1/2" or 1" pebbled Akton® polymer
- reduces shear
- contours easily
- size 16" x 16" to 18" x 18"

### Action® Shear Smart® Pad



- 9/16" cubed Akton® polymer
- reduces surface shearing
- design enhances aeration
- sizes 16" x 16" to 20" x 20"
- custom sizes available to 24" max
- 2 year warranty
- sold standard with Shear Smart® cover



[www.actionproducts.com](http://www.actionproducts.com)

Hagerstown, Maryland 21740 • USA • 301.797.1414

# Order OnLine

Visit [www.actionproducts.com/homecare](http://www.actionproducts.com/homecare) today. Place your order or browse our complete online catalog at your convenience. Here you can:

- read detailed product descriptions
- view all available product options
- download and read clinical studies
- download and read client testimonials
- download and view individual product sheets
- search for a local dealer

Our secure website gladly accepts:



ACTION, AKTON, BODY SMART, SHEAR SMART, XACT and Making a Material Difference are registered trademarks of Action Products, Inc., USA. Centurian, Commuter, Pilot and Professional are trademarks of Action Products, Inc., USA. © 2007



## Action® Hand Exerciser

- 100% Akton® polymer
- improves hand strength
- helps relieve stress
- unlimited degrees of resistance
- easy to heat or cool



## Action® Heel & Elbow Pads

- protects and comforts bony areas
- users can continue daily activities
- users can lie or sleep comfortably
- sold in pairs only
- 100% Akton® polymer



## Action® Commode Pads


- protects and comforts user
- hassle free attachment & care
- closed oval or open split designs
- custom sizes available with diagram



## Action® Comfort Heel Cups

- 100% Akton® polymer
- comforts and supports heels
- very lightweight
- sold in pairs only

[www.actionproducts.com](http://www.actionproducts.com)  
800.228.7763 • 301.797.1414



[www.actionproducts.com](http://www.actionproducts.com)

Action Products, Inc.  
954 Sweeney Drive  
Hagerstown, MD 21740 USA

